

**FUJI COMPUTED RADIOGRAPHY FOR MAMMOGRAPHY (FCR_m)
MAMMOGRAPHY EQUIPMENT EVALUATION (MEE)**

Initial Accreditation _____ Adding New Units to Existing Accreditation _____

FACILITY _____ MAS _____ Model _____ Room _____

MEE Date _____ Within 6 months? _____ PASS/FAIL? _____ Signed by _____

QC Manual Version at facility _____

- 1) Contrast-To-Noise Ratio (CNR) (Tech weekly test)**
 CNR must not vary by more than $\pm 20\%$ of baseline Yes _____ No _____
- 2) Image Quality Evaluation**
Phantom Yes _____ No _____
 4 largest fibers, 3 largest speck groups and 3 largest masses are visible
Phantom Image Quality scores:

Fibers _____	Hard Copy: Background Background OD (>1.20) _____	Soft Copy S value: baseline + 20%
Specks _____	OD ± 0.20 _____	mAs: baseline + 15%
Masses _____	DD ± 0.05 _____	

Printer Quality Control: Follows Manufacturer Specification Yes _____ No _____
Monitor Quality Control: Follows Manufacturer Specification Yes _____ No _____
- 3) S Value Confirmation**
 Corrected S Value should not exceed the range of $120 \pm 20\%$ Yes _____ No _____
- 4) Evaluation of System Resolution** Yes _____ No _____
 Measured performance 8 ± 2 lp/mm, both scan directions
- 5) CR Reader Scanner Performance** Yes _____ No _____
 "T" test image smooth and sharp
- 6) Imaging Plate Fog** Yes _____ No _____
 [is now a semi-annual Tech. test]. Coin not visible on image
- 7) Mammographic Unit Assembly Evaluation** Yes _____ No _____
 Follows Mammographic Unit Assembly Evaluation Form. [The list of items to be checked in the Fuji QC manual should be included as a part of or supplemented to this form].
- 8) Collimation Assessment** Yes _____ No _____
 Deviation between X-ray field and light field is $\leq 2\%$ of SID
 X-ray field does not extend beyond any side of the IR by more than 2% of SID
 and must cover all the IR on the chest wall side
 Chest wall edge of compression paddle doesn't extend beyond IR by more than
 1% SID or appear on image
- 9) AEC System Performance Assessment** Yes _____ No _____
 Measured performance within acceptable limits as stated in the applicable QC manual
 Δ mAs/density step; within 5-15%, COV reprod. ≤ 0.05 , CNR level for 2 & 6 cm within + 15% of the 4 cm value
- 10) System Artifact Evaluation** Yes _____ No _____
 Examine images for artifacts. No objectionable artifacts on image
- 11) Dynamic Range** Yes _____ No _____
 Measured performance within acceptable limits. Exposures with acrylic thicknesses of 0 cm, 2 cm, & 6 cm are discernible
- 12) Primary Erasure** Yes _____ No _____
 Measured performance within acceptable limits. No significant artifacts observed

- 13) **Inter-Plate Consistency** Yes _____ No _____
 Variation of mAs must be within $\pm 10\%$
 Variation of SNR must be within $\pm 15\%$
- 14) **kVp Accuracy and Reproducibility** Yes _____ No _____
 Measured average kVp within $\pm 5\%$ of indicated kVp
 kVp reproducibility coefficient of variation ≤ 0.02
- 15) **Breast Entrance Exposure, AEC Reproducibility and Average Glandular Dose** Yes _____ No _____
 Average Glandular Dose for average breast is ≤ 3 mGy (300 mrad) _____ mrad
 Coefficient of Variation for either R or mAs shall not exceed 0.05
- 16) **Beam Quality Assessment (HVL) Measurement** Yes _____ No _____
 HVL meets minimum specifications in the applicable QC manual. HVL (mmAl) \geq kVp/100 for Mo/Mo & Std. Breast
- 17) **Radiation Output Rate** Yes _____ No _____
 Radiation Output Rate is ≥ 800 mR per second @28 kVp (Mo/Mo) for 3 seconds
- 18) **Viewbox Luminance and Reading Room Conditions** Yes _____ No _____
 Mammographic viewbox is capable of a luminance of at least 3000 nit
 Room illuminance (viewbox surface as seen by observer) is 50 lux or less
 Room illuminance (monitor surface) is: ≤ 20 lux or the limit set by the monitor manufacturer (if less than 20 lux, for soft copy reading)

FUJI COMPUTED RADIOGRAPHY FOR MAMMOGRAPHY (FCR_m)
Annual Survey Report for Reaccreditation
MAMMOGRAPHY EQUIPMENT EVALUATION (MEE)

FACILITY _____ MAS _____ Model _____ Room _____

Survey Date _____ Within 6 months? ____ PASS/FAIL? ____ Signed by _____

QC Manual Version at facility _____

1) Physicist's Review of QC Technologist Tests

Weekly

CNR Test Yes _____ No _____

Phantom Image Test Yes _____ No _____

Printer QC Test Yes _____ No _____

Monitor QC Test Yes _____ No _____

Monthly

Visual Checklist Yes _____ No _____

Quarterly

Repeat Analysis Yes _____ No _____

Semi-annually

Compression Yes _____ No _____

IP Fog Test Yes _____ No _____

2) Image Quality Evaluation

Phantom Yes _____ No _____

4 largest fibers, 3 largest speck groups and 3 largest masses are visible

Phantom Image Quality scores:

Fibers _____

Specks _____

Masses _____

Hard Copy:

Background OD (>1.20) _____

OD \pm 0.20 _____

DD \pm 0.05 _____

Soft Copy:

S value: baseline \pm 20%

mAs: baseline \pm 15%

Printer Quality Control: Follows Manufacturer Specification Yes _____ No _____

Monitor Quality Control: Follows Manufacturer Specification Yes _____ No _____

3) S Value Confirmation

Corrected S Value should not exceed the range of $120 \pm 20\%$ Yes _____ No _____

4) Evaluation of System Resolution

Measured performance 8 ± 2 lp/mm, both scan directions Yes _____ No _____

5) CR Reader Scanner Performance

"T" test image smooth and sharp Yes _____ No _____

6) Mammographic Unit Assembly Evaluation

Follows Mammographic Unit Assembly Evaluation Form. [The list of items to be checked in the Fuji QC manual should be included as a part of or supplemented to this form]. Yes _____ No _____

7) Collimation Assessment

Deviation between X-ray field and light field is $\leq 2\%$ of SID Yes _____ No _____

X-ray field does not extend beyond any side of the IR by more than 2% of SID and must cover all the IR on the chest wall side

Chest wall edge of compression paddle doesn't extend beyond IR by more than 1% SID or appear on image

8) AEC System Performance Assessment

Measured performance within acceptable limits as stated in the applicable QC manual Yes _____ No _____

- Δ mAs/density step; within 5-15%, COV reprod. ≤ 0.05 , CNR level for 2 & 6 cm within + 15% of the 4 cm value
- 9) System Artifact Evaluation** Yes _____ No _____
 Examine images for artifacts. No objectionable artifacts on image
- 10) Dynamic Range** Yes _____ No _____
 Measured performance within acceptable limits. Exposures with acrylic thicknesses of 0 cm, 2 cm, & 6 cm are discernible
- 11) Primary Erasure** Yes _____ No _____
 Measured performance within acceptable limits. No significant artifacts observed
- 12) Inter-Plate Consistency** Yes _____ No _____
 Variation of mAs must be within $\pm 10\%$
 Variation of SNR must be within $\pm 15\%$
- 13) kVp Accuracy and Reproducibility** Yes _____ No _____
 Measured average kVp within $\pm 5\%$ of indicated kVp
 kVp reproducibility coefficient of variation ≤ 0.02
- 14) Breast Entrance Exposure, AEC Reproducibility and Average Glandular Dose** Yes _____ No _____
 Average Glandular Dose for average breast is ≤ 3 mGy (300 mrad) _____ mrad
 Coefficient of Variation for either R or mAs shall not exceed 0.05
- 15) Beam Quality Assessment (HVL) Measurement** Yes _____ No _____
 HVL meets minimum specifications in the applicable QC manual. HVL (mmAl) \geq kVp/100 for Mo/Mo & Std. Breast
- 16) Radiation Output Rate** Yes _____ No _____
 Radiation Output Rate is ≥ 800 mR per second @28 kVp (Mo/Mo) for 3 seconds
- 17) Viewbox Luminance and Reading Room Conditions** Yes _____ No _____
 Mammographic viewbox is capable of a luminance of at least 3000 nit
 Room illuminance (viewbox surface as seen by observer) is 50 lux or less
 Room illuminance (monitor surface) is: ≤ 20 lux or the limit set by the monitor manufacturer (if less than 20 lux, for soft copy reading)